

<u>Medikal Cihazlar Üretim – İthalat – Satış ve Teknik Servis Hizmetleri</u>



Technical Specifications

- * The device works at 220Vac, 50Hz.
- * The device has an adjustable digital timer.
- * At the end of the treatment period, the audible warning starts.
- * At the end of the treatment period, the output is automatically cut off.
- * There are eight independent channels
- * The applied current is monitored on the Light bar
- \ast The frequency of the currents supplied by the device can be set to 2-5-10-50-100-200-250 Hz.
- * The pulse duration of applied current can be adjusted.
- * max output current is 100mA
- * 8 patient cable
- 8 pairs of disposable electrodes 5x5cm
- * The device has two years manufacturing, ten years spare parts warranty.