

Medical Devices Manufacture - Import - Sales and Services



## **MEDELSAN 2 CHANNEL IONTOPHORESIS**

- -The device is suitable for simultaneous Hand-Foot and Armpit applications.
- -Each channel can give 1-20 mA current.
- -Timer can be adjusted from the timer unit.
- -Each channel works independently of each other.
- -The applied current is read from the digital milliAmmeter.
- -There is a patient protection circuit.
- -Device is supplied with administration cups and patient cable.
- -The device has a two-year manufacturing and ten-year spare parts warranty.