



MEDELSAN

M E D I C A L E L E C T R O N I C S

Electronic Eng. Tevfik Soral

Medikal Cihazlar Üretim – İthalat – Satış ve Teknik Servis Hizmetleri



AUDIO-VISUAL EMG BIOFEEDBACK SPECIFICATIONS

- A developed technic of proprioceptive neuromuscular Facilitation.
- APPLICATION FIELDS PSYCHOTHERAPY (Psychophysical Disorders)
- APPLICATION FIELDS URİNER İNCONTİNENS
- PHYSIO-THERAPY
- Nöromuscular Training
- Stress and disorders Hemiplegia and paraplegia
- Connected to stress Disentangling of spasticity
- Stress and migraine Development of voluntary movements
- Asthma All muscular spasms
- Chronic anxiety and fear Cerebral Palsy
- Disorders in learning and attention Spasmodic torticollis
- Stammer treatment All distomic diseases Tensionheadache
- Facial paralysis
- Secures an efficient EMG
- Biofeedback with sound , lightbar and mikrovoltmeter.
- Operates with two 9 V batteries or special adaptor.
- The microvoltmeter Sensitivity can be adjusted to the values of Mikrovolt.
- The entrance impedance is 10 Mohm.
- It possesses two pieces of low filtering and two pieces high Filtering andAdjustable notch filters.
- Sound intensity can be adjusted
- It is movable and 1,5 kg weight.
- Two years quarantee.

Adres: Şakacı Sokak No:64/20 Kozyatağı/Kadıköy/ İstanbul

Tel: 0535 736 36 67 0538 642 89 46 0542 561 67 54 Web: www.medelsan.com

E – Posta: info@medelsan.com